

Introduction

The goals of this training are to develop a better appreciation of equipment controls, scatter radiation in the interventional suite, the radiation output capabilities of the angiography unit, and the effect of various parameters on patient and staff dose.

Section 1 – Principles of fluoroscopic systems

Identify the following:

	Function of this control
Spacer device	
Image receptor	
X-ray tube and housing	
Integrated dose display	
Power on	
Emergency stop	
Integrated radiation shielding	
Collision sensors	

Section 2 – Operation of fluoroscopic systems for interventional purposes

Perform each of the following actions:

- Raise the patient table
- Lower the image receptor
- Collimate the X-ray field
- Change the fluoroscopic pulse rate
- Change the acquisition frame rate
- Store a fluoroscopy loop
- Store a LIH image
- Change the magnification mode

Section 3 – Fluoroscopic exposure outputs and high level control options

Rank the following modes of operation in order of lowest (1) to highest (3) potential radiation output.

- High level fluoroscopy (boost mode) _____
- Acquisition/digital subtraction angiography (DSA) _____
- Normal fluoroscopy _____

Perform each of the following actions:

- Activate normal fluoroscopy
- Activate high level fluoroscopy (boost mode)
- Activate acquisition/DSA

Section 4 – Dose management techniques

Considering the setup after patient transfer, what are some basic actions that should be taken before starting the fluoroscopic procedure? Record the reference air kerma rate ($K_{a,r}$) and the patient entrance air kerma rate (EAKR) in the spaces provided.

- Raise patient table to a comfortable working height
- Lower image receptor as much as practical
- Take one small step down the table or back from the table
- Collimate to the area of interest

	<u>$K_{a,r}$</u>	<u>EAKR</u>
Baseline:		
<u>Order</u>	<u>$K_{a,r}$</u>	<u>EAKR</u>

What impact will the following actions have on the patient EAKR? Will the patient EAKR decrease, increase, or remain the same? Record the $K_{a,r}$ and the EAKR in the spaces provided.

- Action
- Increase magnification
- Reduce the fluoroscopic pulse rate
- Increase in patient thickness

	<u>$K_{a,r}$</u>	<u>EAKR</u>
Baseline:	_____	_____
<u>Impact</u>	<u>$K_{a,r}$</u>	<u>EAKR</u>

Collimate to the area of interest
 Use acquisition/DSA instead of fluoroscopy
 Rotate the C-arm to a 30° oblique projection

What are some other techniques that can be used to reduce the patient dose rate?

Section 5 – Basics of radiation protection

Perform each of the following actions:

- Activate Radiation Release Inhibit
- Activate Positioning without Radiation

What impact will the following actions have on your (the operator’s) dose rate? Will your dose rate decrease, increase, or remain the same? Record the operator dose rate in the spaces provided.

Action

Collimate to the area of interest
 Reduce the fluoroscopic pulse rate
 Raise the patient table
 Lower the image receptor
 Take one small step down or away from the table

	<u>Operator dose rate (mSv/hr)</u>
Baseline:	
<u>Impact</u>	<u>Operator dose rate (mSv/hr)</u>

Where should your personal dosimeter be worn?

Section 6 – Ambient radiation levels and auxiliary shielding

On what side of the patient is the operator dose rate lowest when using lateral (cross-table) fluoroscopy? Test your prediction by making measurements.

How does the ambient radiation level compare at different positions around the fluoroscopy lab?

<u>Position</u>	<u>Relative ambient radiation level</u>	<u>Measured dose rate (mSv/hr)</u>
Operator		
First assistant		
Nurse		
Anesthesiologist		
Outside room door (open)		
Outside room door (closed)		
Control console		

How is the ambient radiation level affected by the use of table drapes?

<u>Position</u>	<u>Predicted % of initial intensity</u>	<u>Measured dose rate (mSv/hr)</u>
Operator		
First assistant		
Nurse		
Anesthesiologist		

How is the ambient radiation level affected by the use of a suspended shield?

<u>Position</u>	<u>Predicted % of initial intensity</u>	<u>Measured dose rate (mSv/hr)</u>
Operator		
First assistant		
Nurse		
Anesthesiologist		

How is the ambient radiation level affected by the use of a sterile protective drape?

<u>Position</u>	<u>Predicted % of initial intensity</u>	<u>Measured dose rate (mSv/hr)</u>
Operator		
First assistant		
Nurse		
Anesthesiologist		

Considering the operator, measure the occupational dose rate with no auxiliary protection and the impact of adding each of the following auxiliary protective devices in sequence.

<u>Protective device(s)</u>	<u>Measured dose rate (mSv/hr)</u>
None	
Table drapes	
Table drapes + suspended shield	
Table drapes + suspended shield + sterile protective drape	

Section 7 – Personal protective equipment

What types of personal protective equipment are available?

What is the impact of a protective garment on occupational dose rate?

What are some considerations for the selection of personal protective equipment?

What characteristics of a protective garment influence its protective value?

What characteristics of protective eyewear influence its protective value?

What is the approximate reduction in occupational dose rate offered by a 0.50 mm protective garment?
A 0.35 mm protective garment?

		<u>mSv/hr</u>
	Baseline:	
<u>Garment</u>		
	<u>Predicted % reduction</u>	<u>mSv/hr</u>
Single 0.35 mm garment		
Single 0.50 mm garment		
Double 0.50 mm garment		

Section 8 – Procedures for recording patient dose data

Where is the air kerma reference point located on your fluoroscope?

Where can you find the following information?

<u>Information</u>	<u>Location</u>
Patient ID	
Type and date of examination	
System ID	
Reference air kerma ($K_{a,r}$)	
Dose or kerma area product (DAP or KAP)	
Fluoroscopy time	
Peak skin dose (PSD)	

Match the dose metric with its description.

<u>Dose metric</u>	<u>Use</u>
$K_{a,r}$ _____	A. Most closely correlated with risk of skin injury.
KAP _____	B. Preferred for notification levels as it is widely available and correlated with skin dose.
Fluoroscopy time _____	C. Best describes total energy imparted to patient, most closely correlated with occupational dose.
PSD _____	D. Useful for quality improvement but poorly correlated with skin dose and occupational dose.

What is a typical value for fluoroscopy time and $K_{a,r}$ for the procedures performed using your fluoroscope? What would be considered exceptional values?

<u>Metric</u>	<u>Typical value</u>	<u>Exceptional value</u>
Reference air kerma ($K_{a,r}$)		
Fluoroscopy time		

At what $K_{a,r}$ should you be concerned regarding the possibility for a skin reaction? Would you expect this value to be reached in your practice?

	<u>Value of concern</u>	<u>Reached in your practice?</u>
Reference air kerma ($K_{a,r}$)		

end