



Fluoroscopic Safety, LLC
sales@fluorosafety.com

Pediatric patient checklist

- Anti-scatter grid removed if patient weighs less than 50 lbs.
- If available, reduced dose pediatric program selected, including pediatric AEC control curve
- Lowest practicable fluoroscopic pulse rate and DSA frame rate selected
- Use of digital acquisition restricted or eliminated, use Store Monitor/Store Fluoro
- Patient as far from X-ray source as practicable
- Image receptor as close to the patient as practicable
- Lowest magnification practicable selected
- X-ray beam collimated tightly to area of interest